

DOCUSTAT DS-220:

- for visualizing traces of writing impressions



Due to the rather large dimensions of the examination unit (450x440 mm), the Docustat DS-220 system can be used for small and large examination objects. The supplied plastic film, type IW, allows the system to be used for a variety of applications such as: Securing writing impression marks (grooves) on paper-like surfaces.

- Simple toner emptying thanks to additional tilting/catching device
 - Faithful reproduction of the original
 - Only reacts to latent traces of lettering - visible lettering is ignored
- Non-destructive - the document is not contaminated
- Simple document preparation with the moistening unit, to increase sensitivity during the electrostatic process
 - Timer for corona development
 - Corona voltage setting
 - Fan - prevents toner particles in the air from flying around

Application example: Print grooves /
Copy marks

A threatening letter was addressed to the US
Trade Center

in Stockholm. A writing pad with the same
form was found in the house of a suspicious
person's house.


With the Docustat system, the print marks could
be made

made visible, and thus the suspect could be
person could be proven to have written the
threatening letter.

could be proven. Depending on how hard the
writing pressure is, the 2nd and 3rd sheets of a
sheets of a bundle of paper or writing pad.
traces can be made visible.

The trace is recorded with the plastic film,
type IW, and by charging the corona unit.

The developing powder is then inserted and
secured with the transparent adhesive film.



Technical data:
Basic equipment

DOCUSTAT DS-220

Electrostatic system

consisting of:

- Basic unit: 450 x 440 x 210 mm with cascade development and integrated fan (prevents toner the toner from flying around), ON/OFF pump, with connection to corona unit
- Vacuum plate: 424 x 424 mm (A3 - interchangeable)
- Suction area inside the vacuum plate vacuum plate: 350 x 385 mm
- Corona unit with ON/OFF switch
- Plastic film for type IW, 1 roll, size 0.4 x 500 m (incl. roll, size 0.4 x 500 m (incl. cutter) cutter)
- Transparent adhesive film, 2 packs 24 pieces each, size 25 x 35 cm
- Cascade developer, 1 kg
- Toner with spoon, 1 kg
- Empty tin

Recommended

Accessories/consumables

Drip tray for

Cascade developer

- Transparent adhesive film, size 25 x 35 cm, 24 pieces (A4)
- Transparent adhesive film, size 50 x 35 cm, 24 pieces (A3)
- Plastic film for imprint marks, type IW, 1 roll, size 0 roll, size 0.4 x 500 m
- Cascade developer, toner



Karl Schröder KG
Mierendorffstrasse 28
D-69469 Weinheim

Web: www.schroeder-prueftechnik.de
@: info@schroeder-prueftechnik.de
T: +49 6201 9068-0
F: +49 6201 9068-29

OPTIONS:

- Document humidification unit
- Digital moisture meter
- Toner pads

Weight: Base unit: 21 kg

Power supply: 100 - 240 Volt

Vacuum pump: -800 mbar

High voltage power supply: 2.5 W



Karl Schröder KG
Mierendorffstrasse 28
D-69469 Weinheim

Web: www.schroeder-prueftechnik.de
@: info@schroeder-prueftechnik.de
T: +49 6201 9068-0
F: +49 6201 9068-29